

## Ortho Clinical Diagnostics

### Current Concepts in HIV Diagnostics

Speaker  
Speaker Title  
Ortho Clinical Diagnostics

This promotional educational activity is brought to you by Ortho Clinical Diagnostics, Inc. and is not certified for continuing medical education. The speaker is compensated by, and presenting on behalf of, Ortho Clinical Diagnostics, Inc. and must present information in accordance with applicable FDA requirements.


\* Product availability subject to local regulatory requirements.

Ortho Clinical Diagnostics

## Agenda

- HIV Infection
- Testing Overview
- CDC's HIV Testing Algorithm
- HIV 4th Generation Combo Tests
- POC HIV Testing
- Summary


## Human Immunodeficiency Virus (HIV)



Source: DOI: 10.13140/RG.2.1.4923.8803. Available from: Jef Hens, Mar 08, 2016

Ortho Clinical Diagnostics

## HIV Life Cycle



Source: DOI: 10.13140/RG.2.1.4923.8803. Available from: Jef Hens, Mar 08, 2016





Ortho Clinical Diagnostics

## Transmission of HIV

- People get or transmit HIV through sexual behaviors and injection drug use
- Mother to child during pregnancy, birth or breastfeeding
- Certain body fluids can transmit HIV:
  - Blood
  - Semen
  - Pre-seminal fluid
  - Rectal fluids
  - Vaginal fluids
  - Breast milk

<https://www.cdc.gov/hiv/basics/transmission.html>

Ortho Clinical Diagnostics



### Global HIV-1 and HIV-2 Epidemiology

- Strains of HIV-1 can be classified into four groups: M, O, N and P
- Group M is responsible for most of the global HIV pandemic
- Group N has only been found in Cameroon
- Group O is responsible for some of the infections in West and Central Africa
- Group P is a new group identified in Cameroon
- 2.7% - 5.4% of the individuals living with HIV worldwide are infected with HIV-2
- Most persons infected with HIV-2 reside in West Africa or immigrate from West Africa
- Fewer than 1% of HIV infections in the United States are caused by HIV-2







Ortho Clinical Diagnostics

### Clinical Manifestations of Acute HIV Infection

- 40-90% develop symptoms of acute HIV<sup>1</sup>
- 50%-90% with symptoms seek medical care
- Of those diagnosed with acute HIV<sup>1</sup>, 50% of patients seen at least 3 times before diagnosis<sup>2</sup>
- The symptoms of HIV that are caused by the body's nature response against the virus include:
  - Fever
  - Fatigue
  - Malaise
  - Myalgia
  - Rash
  - Sore throat
  - Lymphadenopathy<sup>3</sup>

1. Kahn, et al. *New England Journal of Medicine*. 1998:339-33-40.  
2. Weirich A, et al. *Archives of Internal Medicine*. 2003; 163:2097-2100.  
3. Cohen, et al. *Journal of Infectious Diseases*. 2001; 183(1):23-35.


Ortho Clinical Diagnostics



### Comparison of Detection Limits of 4<sup>th</sup> and 5<sup>th</sup> Generation Combination HIV assays

A study supported by NIH, evaluated the sensitivity of HIV clinical diagnostic blood screening and next generation P24 assays:

- Comparison study performed on 20 geographically diverse HIV isolate samples.
- Completed in EQAPOL Laboratories.
- "Similar performance" (sensitivity) across sub-types between 4<sup>th</sup> gen and stand alone p24 Ag assays and Bio-Rad's HIV Ag-Ab assay



- Demonstrated enhanced sensitivity of next generation p24 digital detection platforms, and observed poor performance of POC assays

Stone M, Bainbridge J, Sanchez A et al. Comparison of Detection Limits of 4<sup>th</sup> and 5<sup>th</sup> Generation Combination HIV. Presented at Conference on Retroviruses and Opportunistic Infections, March 2018, Poster #667. Ortho Clinical Diagnostics

### HIV Assay Diagnostic Testing Evolution

Assay progression	1 <sup>st</sup>	2 <sup>nd</sup>	3 <sup>rd</sup>	4 <sup>th</sup>	5 <sup>th</sup>
Year	1985	1987	1991	1997	2015
Antigen (Ag) Source	Virus Infected Cell Lysate	Lysate & Recombinant	Recombinant & Synthetic peptides	Recombinant & Synthetic peptides	Recombinant & Synthetic peptides
Specificity	95-98%	>99%	>99.5%	99.5%	99.5%
Sensitivity	99%	>99.5%	>99.5%	>99.8%	100%
Negative Window	8-10 weeks	4-6 weeks	2-3 weeks	2 weeks	2 weeks
Detects Antibody (Ab) and Ag	IgG Anti HIV-1	IgG anti HIV-1 and IgG anti HIV-2	IgG and IgM anti HIV-1, HIV-2 and Group O	IgG and IgM anti HIV-1, HIV-2 and Group O. Also detects HIV-1 p24 Ag	IgG and IgM anti HIV-1, HIV-2 and Group O. Also detects HIV-1 p24 Ag
Results	Single result	Single result	Single result	Single result; does not differentiate Ab from Ag positivity	Separate HIV-1 and HIV-2 Ab and Ag results
Confirming Tests	HIV-1 western blot (WB) or immunofluorescence (IFA)	HIV-1 WB or IFA, HIV-2 EUSA and WB if HIV-1 confirm is negative	HIV-1 WB or IFA, HIV-2 EUSA and WB if HIV-1 confirm is negative	HIV-1.2 differentiation Assay followed by qualitative HIV-1 RNA PCR if differentiation assay is negative	Not determined at the time of this writing

Alexander, T. S. (2016). Human immunodeficiency virus diagnostic testing: 30 years of evolution. Clinical and Vaccine Immunology 23:249 – 253. doi:10.1128/CVI.00053-16. Ortho Clinical Diagnostics

### Key Dates in the History of HIV Testing

1981	First AIDS case reported
1984	Human immunodeficiency Virus (HIV) identified
1985	First test for HIV licensed (ELISA)
1987	First Western Blot blood test kit
1992	First rapid test
1994	First oral fluid test
1996	First home and urine tests
2002	First rapid test using finger prick
2003	Rapid finger prick test granted CLIA (Clinical Laboratory Improvement Amendments) waiver
2004	First rapid oral fluid test (also granted CLIA waiver)
2006	CDC recommends routine HIV screening in U.S. health care settings
2007	CDC launches Expanded HIV Testing Initiative in U.S.
2007	WHO/UNAIDS global guidelines recommend routine HIV screening in health care settings
2010	First test approved that detects both antigen and antibodies
2012	First rapid oral fluid home test
2012	USPSTF gives routine HIV screening an "A" rating
2013	First rapid test approved that detects both antigen and antibodies, and distinguishes between acute and established HIV-1 infection
2015	Centers for Medicare and Medicaid Services announces Medicare coverage of annual HIV screening for all beneficiaries 15-65, and for those older and younger beneficiaries at "increased risk" for HIV

https://www.kff.org/hiv/aids/topic-sheet/hiv-testing-in-the-united-states/ Ortho Clinical Diagnostics

### General Approach to HIV Diagnostic Testing

- Screening assays
  - Optimized for sensitivity
  - Define samples as negative
- Supplemental/confirmatory assays
  - Optimized for specificity


### HIV Diagnostic Testing Algorithm

HIV Testing is a cornerstone in the fight against HIV spread – recommended in various clinical settings

#### New CDC Recommendations for HIV Testing in Laboratories

*A step-by-step account of the approach*

CDC views recommendations for HIV testing as iterative and applicable to both stand-alone laboratory settings and to point-of-care testing. For more information on the new testing approach, visit <https://www.cdc.gov/hiv/newsroom/2018/04/hiv-testing-recommendations.html>. By putting the latest testing technology to work, laboratories across the United States, we can help address a critical gap in the nation's HIV prevention efforts.



This graphic is designed to illustrate key concepts of the new testing approach to laboratories. For more detail, please see the full guidelines here: <https://www.cdc.gov/hiv/newsroom/2018/04/hiv-testing-recommendations.html>

### FDA-Approved Tests Applicable to the CDC Recommended HIV Testing Algorithm

TEST KIT NAME	MANUFACTURER
<b>Step 1. HIV-1/HIV-2 Ag/Ab combo immunoassay (screening assay)</b>	
Abbott Architect® HIV Ag/Ab Combo Assay	Abbott Laboratories
GS HIV Ag/Ab Combo EIA	Bio-Rad Laboratories
ADVIA Centaur® HIV Ag/Ab Combo	Siemens HealthCare Diagnostics
VITROS® HIV Combo Test	Ortho Clinical Diagnostics
BioPlex® 2200 HIV Ag-Ab Assay	Bio-Rad Laboratories
Roche Elecsys® HIV Combi PT 4th Gen	Roche Diagnostics
<b>Step 2. HIV-1/HIV-2 antibody differentiation immunoassay (supplemental Ab assay)</b>	
Multipot HIV-1/HIV-2 Rapid Immunoassay	Bio-Rad Laboratories
Genieus™ HIV 1/2 Supplemental Assay	Bio-Rad Laboratories
<b>Step 3. HIV nucleic acid test (supplemental RNA assay)</b>	
Apima® HIV-1 RNA Qualitative Assay	Hologic Gen-Probe

https://www.health.ny.gov/diseases/aids/providers/testing/cdc/guidelines.pdf  
 Trademarks not owned by Ortho Clinical Diagnostics are the property of their respective owners. Ortho Clinical Diagnostics

## Guidance for Reporting Results of HIV Diagnostic Testing Algorithm

Step 1 HIV Ag-Ab Screening Assay	Step 2 HIV-1/2 Ab Differentiation Assay	Step 3 HIV-1 RNA Assay	Interpretation for Laboratory Report returned to Health Care Provider	Results Reported
NR	N/A	N/A	Negative for HIV-1 antigen and HIV-1/HIV-2 antibodies. No laboratory evidence of HIV infection.	NOT reportable
R	HIV-1 Pos	N/A	Positive for HIV-1 antibodies. Laboratory evidence of HIV-1 infection is present.	Report all results
R	HIV-2 Pos	N/A	Positive for HIV-2 antibodies. Laboratory evidence of HIV-2 infection is present.	Report all results
R	NR or Any IND (HIV-1, HIV-2 or HIV)	Detected	Positive for HIV-1. Laboratory evidence consistent with acute or early HIV-1 infection is present.	Report all results
R	NR or Any IND (HIV-1, HIV-2 or HIV)	Not Detected	HIV antibodies were not confirmed and HIV-1 RNA was not detected. Possible false positive. Further testing is recommended if warranted by clinical evaluation or risk factors.	Report all results
R	HIV Pos Untypable	N/A	Positive for HIV antibodies. Laboratory evidence of HIV infection is present. Antibodies not differentiated as HIV-1 or HIV-2. HIV-1 RNA and HIV-2 RNA or DNA testing is recommended.	Report all results

<https://www.health.ny.gov/diseases/aids/providers/testing/docs/guidelines.pdf>

Ortho Clinical Diagnostics

## HIV Testing Algorithm Proposed by CDC/APHL

- CDC group tested 228 plasma specimens collected during seroconversion from 26 plasma donors
- Results were analyzed by assay based on current algorithm (WB confirmation) or proposed algorithm (Multispot + NAT confirmation)
- Multispot offered greater sensitivity versus WB
- Combo assay offered higher sensitivity, but only if NAT confirmation was employed

Screening Assay	3 <sup>rd</sup> Generation HIV Antibody Assay	4 <sup>th</sup> Generation HIV Antigen/Antibody Assay
Number of screen positive results	110	135
Western Blot positive	56	56
Multispot positive	90	90
Multispot- or NAT positive	109	134

Trademarks not owned by Ortho Clinical Diagnostics are the property of their respective owners

Source: S Masciotra et al. J Clin Virology 2011; 52S: S17-S22.

Ortho Clinical Diagnostics

## HIV 4<sup>th</sup> Generation Testing

## 4<sup>th</sup> Generation HIV Assay

- Abbott Architect® Ag/Ab Combo
- Bio-Rad BioPlex® Ag/Ab Combo
- Ortho Clinical Diagnostics VITROS® HIV Combo test
- Roche Elecsys® HIV combi PT cobas®
- Siemens Advia Centaur® CHIV

Trademarks not owned by Ortho Clinical Diagnostics are the property of their respective owners

## How to Assess an HIV Test

- Sensitivity: capability of a test to correctly identify people infected with HIV
- Specificity: capability of a test to correctly identify people that are not infected with HIV
- Positive Predictive Value is the probability that a person who tests reactive is indeed infected with HIV
- Negative Predictive Value is the probability that a person who tests negative is not infected with HIV

31


Ortho Clinical Diagnostics

## Antigen Sensitivity

\*\* based on NIBSC/AFSSAPS standards data from three other 4<sup>th</sup> generation manufacturers assay instructions for Use.

**Analytical Sensitivity**  
NIBSC/WHO: 0.46 IU/mL  
AFSSAPS: 12.8 pg/mL

Helps provide assurance in detecting infection



Detection of HIV-1 viral nucleic acid with Nucleic Acid Test (NAT) remains the most sensitive method in identifying acute HIV-1 infection but its use is not widespread due to associated cost, time and labor.

\*VITROS HIV Combo test Instructions for Use, GEM1256\_US\_EN Version 1.0...

32

Ortho Clinical Diagnostics

## VITROS® HIV Combo Test – Clinical Sensitivity\*

**Clinical Sensitivity**  
100% in Adult High Risk

95% Exact Confidence Intervals  
94.94% - 100.00%

**100% in Adult Low Risk**  
95% Exact Confidence Intervals  
91.96% - 100.00%

➔

**Helps Provide Assurance in Detecting Infection**

VITROS HIV Combo test has a clinical sensitivity of 100% in both Adult High and Adult Low Risk Populations

\*VITROS HIV Combo test Instructions for Use, GEM1256\_US\_EN, Version 1.0. Ortho Clinical Diagnostics

## VITROS® HIV Combo Test – Competitive Performance Uncompromised Specificity\*

**Clinical Specificity**

Adult Low Risk 99.58% (CI 99.38%-99.73%)  
High Risk 99.68% (CI 99.06% - 99.93%)

➔

**Fewer false positives can save time, cost and labor in repeat and confirmatory tests**

Testing performed at three external trial sites.  
Assay specificity assessed using 6003 low risk patients and 1004 high risk patients.

High sensitivity while maintaining specificity.

\*VITROS HIV Combo test Instructions for Use, GEM1256\_US\_EN, Version 1.0. Ortho Clinical Diagnostics

## Excellent Precision Can Help Ensure Consistency of Test Results

Panel Member	Mean S/C (Signal to Cut Off)	Within-run CV (%)	Within-lab CV (%)
Anti-HIV-1	0.76	6.6	8.9
Anti-HIV-1	1.28	6.2	7.6
Anti-HIV-1	2.70	5.0	5.8
Anti-HIV-2	0.95	6.6	8.6
Anti-HIV-2	1.30	5.2	7.7
Anti-HIV-2	2.94	5.1	5.8
Anti-HIV-1 Group O	1.14	6.1	7.7
Anti-HIV-1 Group O	1.43	6.8	8.5
Anti-HIV-1 Group O	2.93	5.7	6.6
HIV p24 Ag	0.90	9.1	13.4
HIV p24 Ag	1.58	7.2	10.0
HIV p24 Ag	3.63	6.0	9.2

\*VITROS HIV Combo Instructions for Use, 6842781 and 684 2782, GEM1256\_US\_EN, Version 1.0.  
\*Clinical and Laboratory Standards Institute

Ortho Clinical Diagnostics

## Comparison Between the 4<sup>th</sup> Gen With or Without Results Differentiation

## Seroconversion Panels Data Comparison VITROS® HIV Combo Test vs. BioPlex® 2200

The testing of bleeds from 17 seroconversion panels with both assays shows 100% agreement

PANEL ID	DAYS TO DETECT SEROCONVERSION*		DIFFERENCE
	VITROS HIV Combo test <sup>2</sup> Post bleed day of first reactive result	BioPlex HIV Ag-Ab <sup>1</sup> Post bleed day of first reactive result	
PRB530	0	0	Agreement
PRB569	63	63	Agreement
HIV244	27	27	Agreement
HIV248	18	18	Agreement
HIV9013	23	23	Agreement
HIV9014	0	0	Agreement
HIV9015	30	30	Agreement
HIV9016	30	30	Agreement
HIV9020	89	89	Agreement
HIV9021	47	47	Agreement
HIV9028	53	53	Agreement
HIV9032	22	22	Agreement
HIV9075	22	22	Agreement
HIV9077	42	42	Agreement
HIV9079	40	40	Agreement
HIV9089	16	16	Agreement
HIV12008	28	28	Agreement

\*Seroconversion panels are a group of serial bleeds from plasma donors during seroconversion. They are intended for use by manufacturers and clinical laboratories to evaluate assay sensitivity. Data is from IFU comparison-NOT from a head to head study.  
VITROS® HIV Combo Instructions For Use IFU\_HIV\_Gem1256\_US\_EN  
Bio-Rad BioPlex 2200 System HIV Ag-Ab IFU Ref 685-3455 5007300 July 2015  
Trademarks not owned by Ortho Clinical Diagnostics are the property of their respective owners.

Ortho Clinical Diagnostics

## Florida Department of Health Study

52,000 subjects → 2 cases acute infections confirmed only by nucleic acid testing (Multispot ½ nonreactive, HIV-1 nucleic acid reactive)

FIGURE 2. Reactivity of FDA-approved assays for HIV-1 compared with Western blot.

BM Branson J Acquir Immune Defic Syndr 2010;55:S102-S105. Ortho Clinical Diagnostics

### Clinical Data Comparison VITROS® Systems vs. BioPlex®

- BioPlex 2200 HIV Ag-Ab assay does not shorten window period of HIV infection detection compared to VITROS Systems
- There is no published clinical study that demonstrates the superiority in clinical performance of BioPlex 2200 HIV Ag-Ab assay to 4<sup>th</sup> Gen assays
- Per the current CDC HIV diagnostic testing guidelines, although with separate reporting of antibodies and antigen, BioPlex 2200 HIV Ag-Ab assay is regarded as the 1<sup>st</sup> step assay, and initial reactive samples still require confirmation in the following steps

<https://www.health.ny.gov/diseases/aids/providers/testing/docs/guidelines.pdf>  
Trademarks are the property of their respective owners

Ortho Clinical Diagnostics

### Point-of-Care Testing

### HIV Point-of-Care Testing Facts

- What are HIV Point of Care testing (POCT)?
  - POCTs are screening tests which can be performed on-site while the patient waits and provide results within minutes.
- What is the specimen type for HIV POCT?
  - Oral fluid (saliva), whole blood (finger prick)
- What does an HIV POCT detect?
  - Most HIV POCTs can detect antibodies only; some can detect antibody/antigen
- How long does it take to get a result?
  - The length of time required to get a result will depend on the brand of test used. All HIV POCTs can provide a result in 20 minutes or less

Sources: <https://www.health.ny.gov/diseases/aids/providers/testing/rapid/workbook.htm>

<http://conditions.health.qld.gov.au/HealthCondition/condition/14/116/812/hiv-point-of-care-testing> Ortho Clinical Diagnostics

### HIV Point-of-Care Testing Facts

- What will be the next step if the result is reactive?
  - The reactive result is only preliminary and must be followed-up by confirmatory tests
- What happens if the result is non reactive?
  - No further testing is required. If the patient has a history of recent high-risk activity, repeat testing is recommended

<https://www.health.ny.gov/diseases/aids/providers/testing/rapid/workbook.htm>  
<http://conditions.health.qld.gov.au/HealthCondition/condition/14/116/812/hiv-point-of-care-testing>

Ortho Clinical Diagnostics

### The Clinical Performance of Different HIV POCTs in High Risk Population

- A 3<sup>rd</sup> Gen EIA assay and pooled NAAT were used to provide final results
- Although the performance of different POCTs vary, they are all less sensitive than EIA assay, including the combo rapid assay
- Testing on oral fluids is less accurate than on whole blood
- Currently there is no adequate POC substitute for laboratory-based 4<sup>th</sup> gen testing or HIV NAAT I high-incidence populations

# tests	Sensitivity (95% CI) compared to all cases	Sensitivity (95% CI) compared to EIA+ cases	Specificity (95% CI)
OralQuick® (oral fluid)	2180 51.68 = 75.0% (63.0-84.7)	5160 = 85.0% (73.4-92.9)	2109/2112 = 99.86% (99.59-99.97)
OralQuick (fingerstick)	2175 53.68 = 77.9% (66.2-87.1)	5360 = 88.3% (77.4-95.2)	2107/2107 = 100% (99.82-100)
Uni-Gold™	1614 45.51 = 84.9% (72.4-93.3)	4547 = 95.7% (85.5-99.5)	1561/1561 = 100% (99.76-100)
INSTP™	599 11.15 = 73.3% (44.9-92.2)	1133 = 84.6% (54.6-98.1)	543/544 = 99.82% (98.98-100)
Determine™ Combo	1523 34.40 = 84.6% (70.2-94.3)	3336* = 91.7% (77.5-98.2)	1468/1483 = 98.99% (98.34-99.43)
GS™ HIV-1 HIV-2 Plus O antibody (EIA)	2161 58.66 = 87.9% (77.5-94.6)		2091/2095 = 99.81% (99.51-99.95)

\* This numerator and denominator do not include the participant who tested EIA negative but p24 Ag positive on Determine Combo. Note: estimates cannot be directly compared, as not all POC tests were used on all participants.

<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4762749/pdf/nihms750885.pdf>  
Trademarks not owned by Ortho Clinical Diagnostics are the property of their respective owners

Ortho Clinical Diagnostics

### Considerations for HIV Rapid Tests in New Algorithm

#### Current CDC Recommendation

- FDA approved rapid tests can be used as initial screening test. If any instrumented antigen/antibody test is available, it is preferred due to their superior sensitivity for detecting HIV during acute infection. <https://stacks.cdc.gov/view/cdc/48472>
- For reactive (preliminary positive) result from any rapid HIV test, laboratories should begin confirmation by testing with one of the laboratory-based HIV Ag/Ab Combo immunoassays and follow the steps in the Recommended HIV Laboratory Testing Algorithm

45

Ortho Clinical Diagnostics

## Summary

### Summary

- A fourth generation HIV test detects simultaneously both HIV-1 and HIV-2 antibodies and the HIV p24 antigen
- Provides earlier detection that can help lead to improved diagnosis and prevention
- Current CDC HIV laboratory algorithm recommends the 4<sup>th</sup> Gen HIV test as the initial screening assay
- Although the performance for rapid tests has improved, they cannot replace the initial screening assay in the recommended algorithm

PR-03662

Ortho Clinical Diagnostics

### HIV/AIDs : A manageable chronic disease with advances in treatment

- **Advances in Antiretroviral therapy renders AIDs a chronic and treatable disease more like diabetes.....**

"People living with HIV but benefiting from the latest medical developments can hope to lead normal lives in many respects..."  
The World Health Report, 2003, The World Health Organization

As laboratorians we can play an important role in helping to connect patients to therapy

Source: <http://www.who.int/whr/2003/chapter3/en/>  
<http://www.who.int/publications/10-year-review/hiv/en/index5.html>

Ortho Clinical Diagnostics

## Questions?

## Thank You

© 2018 ORTHO CLINICAL DIAGNOSTICS  
ALL RIGHTS RESERVED

PR-03662