



The Art of Quality Control, Selected Topics

Laboratory Education Expo, April 2019

Michael Toyoshima, BSMT (ASCP)SC, CLS

mike_toyoshima@bio-rad.com 1-800-854-6737, x1272



Objectives

- How to reliably derive QC statistics
- Use resources to determine correct data configuration
- Explain selected Westgard Rules and their behavior
- Consider different ways to capture blood bank, microbiology and infectious disease QC

2



Most FAQ...Why Does QC Matter?

Point-of-care testing for HIV: HIV counselling and testing

[BJ Johnston](#), MD FRCP¹ and [JM Conly](#), MD CCFP FRCP FACP²

HIV Testing After a First Positive Rapid Diagnostic Test: A Role for Nucleic Acid Testing? 

[Anne M Neilan](#) , [Jennifer E Cohn](#), [Jean-Francois Lemaire](#), [Emma Sacks](#), [Rebecca Alban](#), [Kenneth A Freedberg](#), [Rochelle P Walensky](#), [Andrea L Ciaranello](#)

Open Forum Infectious Diseases, Volume 5, Issue 8, August 2018, ofy170,
<https://doi.org/10.1093/ofid/ofy170>

Published: 29 August 2018 [Article history](#) ▼

3



Arriving at Decision Limits - New Control Testing

- Account for Variability
- Detect and Eliminate Outliers
- Calculate a +/- 2 SD Range (~ 95% CI*)

*CI = Confidence Interval

BIO-RAD

Initial Control Value Testing

Philosophy: Have as much variety in test conditions as seen in daily testing

- Personnel
- Testing Parameters
 - Reagents
 - Calibrations
 - Environment (Temperature & Humidity)
 - Electrical Supply


BIO-RAD

Gathering Data Points

The figure shows a normal distribution curve with the following data points:

Standard Deviation Interval	Percentage of Occurrence
Between -1s and +1s	68.0%
Between -2s and +2s	95.5%
Between -3s and +3s	99.7%

BIO-RAD



Coefficient of Variation [CV], Uses

Compare methods, reagents and instruments under similar conditions

Analyte Sample Level of Analyte

Use to Compare:

- Proficiency Survey Results
- Manufacturer's Claims
- Monthly QC Peer Group Reports
- Your Own Historical Performance

BIO-RAD

Effect of Mean Value on CV

HDL (mg/dL)			
n	80	80	80
Mean	20.6	51.4	128.7
Within run SD	1.1*	← SD Similar →	1.4*
Within run CV (%)	5.5	2.6	1.1
Total SD	1.4	1.4	2.8
Total CV (%)	6.9	2.8	2.2

* SDs similar

BIO-RAD

Effect of Mean Value on CV

HDL (mg/dL)			
n	80	80	80
Mean	20.6	51.4	128.7
Within run SD	1.1	1.3	1.4
Within run CV (%)	5.5*	2.5	1.1*
Total SD	1.4	1.4	2.8
Total CV (%)	6.9	2.8	2.2

* CVs 5X different due to different means

10

Coefficient of Variation Ratio [CVR]*

Calculation

$$\frac{\text{Lab's Monthly CV}}{\text{Peer Monthly CV}}$$

Ideally, CVR ≤ 1.0, since your values are from a single lab, while the peer CV is from several.

If CVR = 1.5 to 2.0, the lab is 50-100% less precise than its peer group, usually requiring investigation.

*Also known as CVI, Coefficient of Variation Index

11

Coefficient of Variation Ratio [CVR]

Calculation

$$\text{CVR} = \frac{5.0}{5.0} = 1.0$$

Your imprecision is same as peers

$$\text{CVR} = \frac{10.0}{5.0} = 2.0$$

- What is the probability? 1 out of ? Labs?
- Can the CVR be too 'small'? Answer later

12

Standard Deviation Index, SDI

- Calculation:
- SDI = $\frac{[\text{Lab Mean} - \text{Peer Group Mean}]}{\text{Peer Group's 1 standard deviation}}$
- Use to assess bias compared to peer group

TARGET SDI = 0.0, so lab's mean value is the same as the peer's \rightarrow NO bias

13 BIO-RAD

Standard Deviation Index, SDI

CALCULATION

SDI = $\frac{100 - 100}{1.0} = 0.0$

Your mean has no bias against peers

SDI = $\frac{120 - 100}{10.0} = 2.0$


What is the probability? 1 out of ? Labs?
Can the SDI be too small? Answer later.


14 BIO-RAD


Textbook Example #1 – Imprecision versus Bias

GOOD PRECISION
NO BIAS

15 BIO-RAD







Imprecision Present – Bias Present?

For troubleshooting, verify precision first

19 BIO-RAD

Traditional Calculation of +/- 2SD Range

Calculation of A 95% CI

Estimation of a 2s Range for a small sample by calculation a 95% confidence interval.

Data Set: Level II — Fibrinogen

408	397	400	400	
402	393	391	389	
401	403	393	385	n = 30
401	404	393	389	\bar{x} = 396.7
394	404	389	387	s = 7.94
402	418	385	392	2s range = 380.8 – 412.6
392	400	411		
392	388	398		

20 BIO-RAD

Traditional Calculation of +/- 2SD Range

Calculation of A 95% CI

Estimation of a 2s Range for a small sample by calculation a 95% confidence interval.

Data Set: Level II — Fibrinogen

408	397	400	400	$CV = 7.94/396.7 \times 100\% = 2.0\%$
402	393	391	389	
401	403	393	385	n = 30
401	404	393	389	\bar{x} = 396.7
394	404	389	387	s = 7.94
402	418	385	392	2s range = 380.8 – 412.6
392	400	411		
392	388	398		

**IF +/- 2 SD, total SD = 4. Divide by 4 to find 1 SD.
(1 SD/mean)/100% = CV%. Compare with Insert Range**

21 BIO-RAD

In-Lab Resources, Manuf Controls

- Manufacturer Controls, Bio-Rad Mfr Reports

Roche/Hitachi cobas c 502 analyzer

Value sheet

Short name / component	Methods	ACN	Value	Range	1s	Unit
NH3L Ammonia	enzymatic	8478	227	182 - 272	15	µmol/L
	CV = (15/227) x 100% = 6.6%					

Unity

Manufacturer Report for Roche

Ethanol/Ammonia • Lot 54270 • Exp 31-Oct-2020

Roche cobas 6000/8000c 311	Level	Mon	Cum	Level	Mon	Cum
Mean	41.12	39.93	110.2	105.3	285.4	289.8
SD	3.26	3.22	4.8	4.86	5.12	5.26
CV	7.7	8.1	4.3	4.6	3.2	3.3
# Points	2627	4984	1185	2562	1727	2912
# Labs	85	85	38	38	85	75

CV ~ 3-5% at comparable concentration

25

BIO-RAD

Look for Imprecision Claims

Precision Roche NH3L, 'c' models

Precision was determined using human samples and controls in an internal protocol with repeatability* (n = 21) and intermediate precision** (3 aliquots per run, 1 run per day, 21 days). The following results were obtained:

Repeatability*	Mean	SD	CV
	µmol/L (µg/dL)	µmol/L (µg/dL)	%
AEC Control N	60.7 (103)	1.4 (2)	2.3
AEC Control A	202 (344)	2 (3)	0.8
Human plasma 1	28.6 (48.7)	2.5 (4.3)	8.8
Human plasma 2	585 (996)	1 (2)	0.2
Intermediate precision**	Mean	SD	CV
	µmol/L (µg/dL)	µmol/L (µg/dL)	%
AEC Control N	56.9 (97.1)	2.8 (4.8)	4.9
AEC Control A	203 (346)	4 (7)	1.8
AEC Control N 1.2 dil.	28.1 (47.7)	2.6 (4.4)	9.4
AEC Calibrator	318 (542)	5 (9)	1.5

* repeatability = within-run precision

** intermediate precision = total precision / between run precision / between day precision

26

BIO-RAD

Group Values by Method

Unity

Worldwide Report

Unassayed Chemistry • Lot 31880 • Exp 31-May-2020

Glucose	Hexokinase	mg/dL	Level	Mon	Cum	Level	Mon	Cum
Beckman AU	Beckman Coulter AU Series	85 / 284	1	85.03	85.08	2	283.8	283.8
	Mean			1.82	1.85		5.75	6.37
	SD			2.1	2.2		2.0	2.2
	CV			5333	25491		5327	25800
	# Points			137	154		138	155
	# Labs							
Roche cobas	Roche cobas 6000/8000c 311	85 / 283	1	84.63	84.78	2	282.4	282.9
	Mean			1.61	1.70		5.08	5.30
	SD			1.9	2.0		1.8	1.9
	CV			16945	84266		16667	83786
	# Points			336	389		334	387
	# Labs							
Hexokinase	Method Group • Hexokinase	85 / 282	1	84.46	84.71	2	281.6	282.4
	Mean			2.25	2.18		6.13	6.23
	SD			2.6	2.6		2.2	2.2
	CV			38294	2039K		38941	2039K
	# Points			913	1044		913	1042
	# Labs							

Group Values by Method

27

BIO-RAD

Group Values by Method, Hematology

Unity Worldwide Report
Hematology-16 • Lot 77800 • Exp 28-Apr-2019

Models from Abbott to Sysmex

Electrical Impedance	Level	Mon	Cum	Level	Mon	Cum	Level	Mon	Cum
Method Group • Electrical Impedance									
Mean	1	2.29	2.32	2	4.21	4.20	3	5.15	5.15
SD		0.081	0.101		0.113	0.129		0.156	0.160
CV		3.5	4.4		2.7	3.1		3.0	3.1
# Points		1914	4120		2381	6794		2322	6610
# Labs		136	149		160	224		153	211

Hundreds of labs, different manufacturers/models, and **%CV < 5%**
ALL using method ELECTRICAL IMPEDANCE

28 BIO-RAD

Your historical performance

Unity Laboratory Histogram
Multiquel 1.2.3 • Lot 45780 • Exp 31-Oct-2019
February 2019 • Lab COBAS 6000 #1

Your historical performance to use for reference

29 BIO-RAD

Use historical CV to set new ranges

Coefficient of Variation, CV

Measure of variability (random error or imprecision)
Expressed as a % value

CALCULATION: $\frac{1 \text{ standard deviation, (s)}}{\text{Mean Value}} \times 100\%$

Example: Fibrinogen $\frac{28 \text{ mg/dL}}{280 \text{ mg/dL}} \times 100\% = 10.0\% \text{ CV}$

CV REMAINS the same; **SD CHANGES** to match new mean; so, now 28 mg/dL is a reasonable 1 SD value

30 BIO-RAD

Good initial control testing – WHY?

Establishing Data Parameters


Determination of valid mean and standard deviation values are crucial to successful data acceptance and data rejection by error detection schemes

31

BIO-RAD

Range Defined After Initial Testing

Data Distribution




32

BIO-RAD

Shift (down) in Values – Systematic Error

Data Distribution

Unexpected Drop in Values



33

BIO-RAD

Erratic Values – Random Error

Data Distribution

Unexpected Imprecision

34 BIO-RAD

Same Mean, 1 SD and CV: Different Stories

'normal'

increase ~ 11th

cyclical

In a month of results, do YOU see any values > 2SD?

35 BIO-RAD

Gaussian Distribution

$\pm 1S - 68.0\%$
 $\pm 2S - 95.5\%$
 $\pm 3S - 99.7\%$


1. Probability
 2. Size of Error
 3. Important?

Probability
 $> \pm 1s: 1/3$
 $> \pm 2s: 1/20$
 $> \pm 3s: 1/333$

36 BIO-RAD

1_{2s} RULE


- Warning rule only: look for errors in the current run and in previous runs
- Any control value greater than two standard deviations from the mean
- Does not require run rejection



1_{2s} RULE

- Warning rule only: look for errors in the current run and in previous runs
- Any control value greater than two standard deviations from the mean
- Does not require run rejection


● Level 1 ■ Level 2



• Level 1 & 2 LJ's superimposed


1_{2s} RULE

- Warning rule only: look for errors in the current run and in previous runs
- Any control value greater than two standard deviations from the mean
- Does not require run rejection



1_{3s} RULE


- Any control value greater than three standard deviations on either side of the mean
- Detects principally random error
- May indicate the beginnings of large systematic error
- Applied in a single run
- May be a cause for run rejection



- 1 of 2 Westgard rules for imprecision


2_{2s} RULE (Across Run)

- Two Consecutive Control values each from separate runs between two and three standard deviations from the mean both on the same side of the mean.
- Detects systematic error which is potentially affecting a specific portion of the method curve.
- May be a cause for run rejection



2_{2s} RULE (Within Run)

- Two control values between 2 and 3 standard deviations from the mean, both on the same side of the mean.
- Detects systematic error which is potentially affecting the entire method curve.
- May be a cause for run rejection




- Level 1 & 2 LJ's superimposed

Any labs using this Westgard Rule?

10 \bar{x} RULE (Within Control Materials)

- Ten consecutive values for the same level of control on the same side of the mean.
- Indicates a bias in the test system potentially affecting a specific portion of the method curve.
- Usually is not a clinically significant bias.
- Run rejection is not required.



Frequently Asked Questions

My control values do not fall within your package insert range. What should I do?

- Obtain peer information online
- Any changes to instrument, reagents/calibrators, software? Any manufacturer notifications?
- Call your QC program
- Confirm change affects patients / QC? Values?

44

BIO-RAD

Bad Habits of Quality Control

- Repeat the Control
 - If **out** 1 out of 20, then **in** 19 out of 20, right?
 - If rules chosen properly, less need to repeat
 - See next section
- Open a New Bottle of Control
 - Proper preparation & storage – training
 - Expensive
- Recalibrate
 - Introduces bias, may mask other problems

From "QC - The Out of Control Problem", Elsa Quam, Westgard.com

45

BIO-RAD

A New View of 1-2s Rule Usage

• Repeat the Control
 – If out 1 out of 20, then in 19 out of 20, right?

Clinical Chemistry 58:5
925-929 (2012)

Should I Repeat My 1:2s QC Rejection?
John A. Korte¹, Lubert Lachant², Leslie L. Cook-Richard³

Introduction: Repeating QC that is outside 2SD from the mean is the only alternative to the common practice of discarding the control material and replacing it with a new lot. This approach has been widely debated.


Options: Two sampling strategies are discussed: (1) repeat the control material and (2) discard the control material and replace it with a new lot. The latter strategy is generally preferred because of the lower risk of a false rejection.

Conclusion: The decision to repeat or discard the control material should be based on the risk of a false rejection and the cost of a false rejection. The latter strategy is generally preferred because of the lower risk of a false rejection.

AACC:
 Clinical Chemistry
 58:5; pg 925-9
 2012

<http://clinchem.aaccjnl.org/content/clinchem/58/5/925.full.pdf>

46 BIO-RAD




Qualitative Responses Selected Considerations

- **Urinalysis**
 - Added by Analyte/Instrument/Strip
 - In the past, Reagent = Dedicated
 - Now, designate dipstick 'name'
 - 10SG, 9SG = same responses
 - May differ when manual vs instrument
 - If so, need to create new combo
 - May differ geographically
 - **Would like to compare, if possible**

48 BIO-RAD


Qualitative Responses Selected Considerations

- **Microbiology**
 - Antimicrobial Susceptibility Testing (AST)
 - Analyte currently ANTIBIOTIC (µg)
 - Method currently agar
 - Considerations for additives, atmosphere, time, and temperature
 - Future to build MIC (Min Inhib Conc), same considerations, but liquid media
 - Etest MIC, various incubations, reporting in ranges of concentrations
 - Specs may differ by source/site of specimen

49 


Qualitative Responses Selected Considerations

- **Microbiology, continued.**
 - AST testing scheme, 'Lot' & AMIRT
 - **LOT, selected as ATCC (American Type Culture Collection)**
 - Analyte, the antibiotic (AST or MIC)
 - Method, the growth conditions
 - Instrument, manual or automated, if significant
 - Reagent, disk or broth manufacturer
 - Temperature, growth condition
 - Unit, **quantitative (mm)**

50 

Qualitative Responses Selected Considerations

- **ImmunoHematology (IH = Blood Bank) QC**
 - Covers 3 common tasks: ABO-Rh; Ab Screening/Ab Identification
 - Even within ABO-Rh, there are multiple combinations of Ab-Ag testing, especially the Rh subgroups; QC testing of BB materials is done daily
 - Initially, the QC will be a documentation of daily QC testing
 - Eventually, hope to have applications for troubleshooting marginal reagents and evaluating changes in titer for Ab and Ag

51 

Qualitative Responses Selected Considerations

- **Infectious Disease Testing**
- These tests are normally reported qualitatively as Neg/Nonreactive, Pos/Reactive; sometimes Indeterminate.
- Reporting standards – per manufacturer – are usually based on Sample/Cutoff ratios, which are measured quantitatively
 - 'Reactive' may be based on S/CO of > 1.00
 - 'Reactive' may be based on S/CO of < 1.00

52

BIO-RAD

Infectious Disease Testing

- How to monitor QC of testing that is reported as Reactive or Non-Reactive?
- Look at the insert:
- First you have RLU's and S/CO:

RESULTS

Calculations

- The ARCHITECT / System calculates the cutoff RLU from the mean RLU of the three Calibrator 1 replicates and stores the result. The cutoff RLU is determined by multiplying the Calibrator 1 Mean RLU by 0.40.
Cutoff RLU = Calibrator 1 Mean RLU x 0.40
- The ARCHITECT / System calculates a result based on the ratio of sample RLU to the cutoff RLU (S/CO) for each specimen and control.

$$S/CO = \text{Sample RLU} / \text{Cutoff RLU}$$

53

BIO-RAD

Infectious Disease Testing


- How to monitor QC of testing that is reported as Reactive or Non-Reactive?
- Look at the insert:
- Consider the interpretation

Interpretation of Results

ARCHITECT HIV Ag/Ab Combo Initial Result		
Initial Result (S/CO)	Instrument Interpretation	Retest Procedure
< 1.00	NONREACTIVE (NR)	No retest required.
≥ 1.00	REACTIVE (R)	Retest in duplicate.

54

BIO-RAD



The Art of QC, Selected Topics

Questions?

56 BIO-RAD
