

Genetics 101 for BRCA and Non-Invasive Prenatal Testing:

Presented by:
Steven Keiles, MS, LCGC
Director, Genetic Counselor Organization

September 28, 2015



Presentation Agenda

- Cancer Genetics 101
- "Red Flags" for hereditary cancer
- *BRCA1* and *BRCA2* genes
 - Cancer risks, testing criteria, and management guidelines
- Additional hereditary breast cancer genes
 - Cancer risks, testing criteria, and management guidelines
- Noninvasive prenatal screening (NIPS)
 - Chromosomal abnormalities
 - Evolution of prenatal screening
 - Clinical utility of NIPS
- Genomic Client Services

CONFIDENTIAL: Not for distribution, copyright reserved Quest Diagnostics 2015



Cancer Genetics 101

CONFIDENTIAL: Not for distribution, copyright reserved Quest Diagnostics 2015



How Does Cancer Develop?


- Sporadic cancer
 - Sporadic cancer develops due to genetic changes within the cells of a specific body tissue
 - Causes for sporadic cancer include aging, hormone history, radiation exposure, smoking, alcohol use, BMI
 - **Most contributing factors are still unknown**
- Hereditary cancer
 - Hereditary cancer results from a mutation in a gene
 - Gene mutations can be passed down in families, causing an increased susceptibility to cancer
- All cancer is **genetic**, but most cancers are not **hereditary**

CONFIDENTIAL: Not for distribution, copyright reserved Quest Diagnostics 2015



How Does Cancer Develop?

Hereditary



Sporadic



CONFIDENTIAL: Not for distribution - copyright reserved Quest Diagnostics 2015



Examples of Genetic Variation

Color → Colour

- ❖ changed spelling but not meaning of word
- ❖ benign polymorphism


Hear → Ear

- ❖ letter deletion changes the word
- ❖ pathogenic mutation

CONFIDENTIAL: Not for distribution - copyright reserved Quest Diagnostics 2015



What Percentage of Cancer is Hereditary?



CONFIDENTIAL: Not for distribution - copyright reserved Quest Diagnostics 2015



“Red Flags” for Inherited Susceptibility to Cancer


- Cancer in 2 or more closely related relatives
- Multiple generations affected
- Early age at diagnosis
- Multiple primary tumors
- Bilateral or rare cancers
- Constellation of tumors consistent with a specific cancer syndrome
- Certain ethnic backgrounds (e.g. Ashkenazi Jewish ancestry)



CONFIDENTIAL: Not for distribution - copyright reserved Quest Diagnostics 2015




Assessing Your Patient's Family History



CONFIDENTIAL: Not for distribution - copyright reserved Quest Diagnostics 2015



Autosomal Dominant Inheritance



CONFIDENTIAL: Not for distribution - copyright reserved Quest Diagnostics 2015



Genetic Testing Strategy

- **The right person**
 - When possible, testing should begin with the individual most likely to test positive (family member who has actually had cancer, youngest age at diagnosis, bilateral disease, multiple primary cancers)
- **The right test**
 - Consider:
 - Is there a known familial mutation?
 - Is the patient of Ashkenazi Jewish ancestry?
 - What types of cancers have been diagnosed in the patient and/or the family?
 - Single gene(s) vs targeted panel vs large panel
- **The right time**
 - Are the test results needed for surgical decision making?
 - Consider the psychosocial implications of genetic testing
 - Testing is generally not warranted for individuals under the age of 18

CONFIDENTIAL: Not for distribution - copyright reserved Quest Diagnostics 2015



Hereditary Breast and Ovarian Cancer Syndrome
BRCA1 and BRCA2

CONFIDENTIAL: Not for distribution - copyright reserved Quest Diagnostics 2015

BRCA1 and BRCA2

The most common high-risk hereditary breast cancer genes

- Autosomal dominant
- Carrier frequency ~1/400 in general population
- Carrier frequency ~1/40 in Ashkenazi Jewish population

Cancer Risks to age 70

	BRCA1	BRCA2
Female Breast Cancer	55% to 65%	45% to 47%
Ovarian Cancer	39%	11% to 17%
Male Breast Cancer	1.2%	6.8%

*Also, increased risks for prostate cancer, melanoma, and pancreatic cancer

CONFIDENTIAL: Not for distribution, copyright reserved Quest Diagnostics 2015



Medical Management Guidelines: BRCA positive

Cancer	Intervention	Age and Frequency
Breast	Screening	Clinical breast exam every 6-12 months, starting at age 25 y for women and 35 y for men Annual breast MRI, starting at age 25 y for women, add annual mammogram at age 30; consider baseline mammogram at 40 y for men
	Chemoprevention	(i.e. Tamoxifen, Raloxifene, aromatase inhibitors)
	Surgery	Discuss risk-reducing mastectomy
Ovarian	Surgery	Risk-reducing bilateral salpingo-oophorectomy, ideally between the ages of 35 and 40 y
	Chemoprevention	(i.e. oral contraceptives)
	Screening	Transvaginal ultrasound and CA-125 blood test; repeat every 6 months starting at age 30-35 y*
Prostate	Screening	Starting at age 40 y for BRCA2 carriers / consider for BRCA1 carriers
Melanoma	Screening	No specific guidelines: Consider annual dermatological examinations
Pancreatic	Screening	No specific guidelines: Consider pancreatic cancer screening if + family history

*Ovarian cancer screening has not been shown to be sufficiently sensitive or specific to support routine screening recommendation but can be considered at the clinician's discretion

CONFIDENTIAL: Not for distribution, copyright reserved Quest Diagnostics 2015



Testing Criteria for BRCA1 and BRCA2

Family History of Breast Cancer

- Relative with a previously identified BRCA1 or BRCA2 mutation
- 1st/2nd-degree blood relative who meets any criteria in the Personal History sections
- 3rd-degree relative with breast^a and/or ovarian^b cancer and ≥2 close blood relatives^c with breast and/or ovarian^b cancer

Personal History of Other (Nonbreast) Cancers

- Epithelial ovarian cancer^a
- Pancreatic or prostate cancer with ≥2^d close blood relatives^e diagnosed with breast, ovarian^b, pancreatic, or prostate cancer (Gleason score ≥7)

Personal History of Breast Cancer^a

Age at Diagnosis	Additional Criteria Only 1 of the following is necessary.
≤45 y	<ul style="list-style-type: none"> • No additional criteria necessary
≤50 y	<ul style="list-style-type: none"> • ≥2 primary breast tumors^a • ≥1 close blood relative^e with breast cancer • Limited family history
≤60 y	<ul style="list-style-type: none"> • Breast cancer that is negative for ER, PR, and HER2 (triple negative)
Any age	<ul style="list-style-type: none"> • Patient is male • ≥1 close blood relative^e with breast cancer diagnosed by age 50 or with epithelial ovarian^b cancer diagnosed at any age • ≥2 close blood relatives^e with breast cancer • ≥2 close blood relatives^e with prostate cancer (Gleason score ≥7) or pancreatic cancer • ≥1 close male blood relative^e with breast cancer • Ethnicity (eg, Ashkenazi Jewish) associated with higher mutation frequency

CONFIDENTIAL: Not for distribution, copyright reserved Quest Diagnostics 2015



Genetic Testing Strategy – BRCA1/BRCA2

BRCAvantage™, Comprehensive (includes large rearrangements)

BRCA1

BRCA2

BRCAvantage™, Ashkenazi Jewish Screen

BRCA1

BRCA2

BRCAvantage™, Single Site *

* Except for patients of Ashkenazi Jewish ancestry


CONFIDENTIAL: Not for distribution, copyright reserved Quest Diagnostics 2015



Beyond BRCA:

Additional Breast Cancer Susceptibility Genes


CONFIDENTIAL: Not for distribution, copyright reserved Quest Diagnostics 2015



Beyond BRCA: Other genes

Mutations in *BRCA1* and *BRCA2* explain 15%-20% of hereditary breast cancer cases

Mutations in *TP53*, *PTEN*, *STK11*, *CDH1*, and *PALB2* explain an additional 3% – 4.5% of hereditary breast cancer



CONFIDENTIAL: Not for distribution, copyright reserved Quest Diagnostics 2015

Couch et al., Science 343 (6178): 1466-1470



Beyond BRCA: Gene Panels

BRC Advantage Plus™

BRCA1, BRCA2, TP53, PTEN, STK11, CDH1, PALB2

- Simultaneous analysis of seven genes associated with hereditary breast and other cancers
- These genes all have NCCN guidelines for medical management if results are positive
- More efficient and cost-effective if personal/family history raises concern for more than one syndrome
- May be useful for those who have previously tested negative for a single syndrome, but whose personal or family history remains strongly suggestive of an inherited susceptibility
- Also available as a reflex test if *BRCA* analysis is negative

CONFIDENTIAL: Not for distribution, copyright reserved Quest Diagnostics 2015




BRC Advantage Plus™

Associated Cancer	Gene						
	<i>BRCA1</i>	<i>BRCA2</i>	<i>TP53</i>	<i>PTEN</i>	<i>CDH1</i>	<i>STK11</i>	<i>PALB2</i>
Breast	•	•	•	•	•	•	•
Colorectal			•	•	•	•	
Endometrial			•	•	•	•	
Stomach					•		
Melanoma	•	•	•	•			
Ovarian	•	•	•			•	
Pancreatic	•	•	•				•
Prostate	•	•					
Other			•	•		•	

CONFIDENTIAL: Not for distribution, copyright reserved Quest Diagnostics 2015



Breast Cancer Risks



* Lifetime risk for cancer of any type, including breast



CONFIDENTIAL: Not for distribution


Genetic Testing Criteria and Management Guidelines

Gene	Genetic Testing Criteria	Management Guidelines
BRCA1/BRCA2	NCCN: Genetic/Familial High-Risk Assessment: Breast and Ovarian	
TP53		
PTEN		
CDH1		
CDH1	NCCN Clinical Practice Guidelines in Oncology: Gastric Cancer	
STK11	None (Clinical Diagnosis)	NCCN: Genetic/Familial High-Risk Assessment: Colorectal
PALB2	None	ACS recommends screening with MRI for women with at least 20% lifetime risk of breast cancer



CONFIDENTIAL: Not for distribution - copyright reserved Quest Diagnostics 2015

BRCAVantage Plus™ Technology



CONFIDENTIAL: Not for distribution - copyright reserved Quest Diagnostics 2015

Interpreting Results of Genetic Tests

- Interpretation Summary
- Color-Coded 5-Tier Classification
 - Known Pathogenic (**RED**)
 - Likely Pathogenic (**RED**)
 - VUS (**YELLOW**)
 - Likely Benign (not reported)
 - No Mutation Detected (**GREEN**)
- Comprehensive Interpretation
- ACMG Guidelines

CONFIDENTIAL: Not for distribution - copyright reserved Quest Diagnostics 2015

Possible Test Results

POSITIVE (Pathogenic, Likely Pathogenic):

- Refer to NCCN Management Guidelines
- Relatives should consider single-site testing

NEGATIVE:

- Uninformative result versus true negative. True negative only if family history is due to a known mutation.
- Consider additional testing for your patient and/or other family members

VARIANT OF UNCERTAIN SIGNIFICANCE:

- Additional data required for interpretation (e.g., family studies, functional studies)
- Patient recommended to check in with MD/Quest Diagnostics on regular basis (e.g., every year) as variant status may change
- Management based on patient's personal and family history

CONFIDENTIAL: Not for distribution. copyright reserved Quest Diagnostics 2015



- Co-founded by Quest Diagnostics and Insemr
- Provides scientists and commercial laboratories around the world open access to *BRCA1* and *BRCA2* genetic data
- The goal of BRCA Share is to accelerate research on *BRCA* gene mutations and clarify the pathogenicity of variants of uncertain significance

CONFIDENTIAL: Not for distribution. copyright reserved Quest Diagnostics 2015



Noninvasive prenatal screening (NIPS)


QNatal™ Advanced



- Overview
 - Quest Diagnostics
 - Chromosomal abnormalities
 - Evolution of prenatal screening
 - Clinical utility of NIPS
- Performance of QNatal Advanced
- Summary

CONFIDENTIAL: Not for distribution. copyright reserved Quest Diagnostics 2015








Chromosomal Aneuploidy

- Common autosome
 - Trisomy 21 (Down syndrome): 47, XX,+21
 - Trisomy 18 (Edwards syndrome): 47, XX,+18
 - Trisomy 13 (Patau syndrome): 47, XX,+13
- Common sex chromosome
 - Monosomy X (Turner syndrome): 45, X
 - Triple X (Triple X syndrome): 47, XXX
 - Klinefelter (Klinefelter syndrome): 47, XXY
 - XYY (XYY syndrome): 47, XYY



Chromosomal Aneuploidy: Clinical Picture

- Common autosome
 - Trisomy 21 (Down syndrome): 47, XX,+21
 - Trisomy 18 (Edwards syndrome): 47, XX,+18



CONFIDENTIAL: Not for distribution, copyright reserved Quest Diagnostics 2015



Chromosomal Microdeletions

- Microdeletion: a chromosomal deletion too small to be detected by standard cytogenetic analysis
- Microdeletion syndrome - Characteristics vary based on the size and location of the deletion
 - DiGeorge syndrome: 22q
 - 1p36 syndrome: 1p
 - Angelman syndrome: 15q
 - Prader-Willi syndrome: 15q
 - Cri-du-chat syndrome: 5p
 - Wolf-Hirschhorn syndrome: 4p
 - Jacobsen syndrome: 11q
 - Langer-Giedion syndrome: 8q



CONFIDENTIAL: Not for distribution, copyright reserved Quest Diagnostics 2015

Chromosomal Microdeletions: Clinical Picture

- Microdeletion: a chromosomal deletion too small to be detected by standard cytogenetic analysis
- Microdeletion syndrome - Characteristics vary based on the size and location of the deletion
 - DiGeorge syndrome: 22q
 - Cri-du chat syndrome: 5p
 - Angelman syndrome: 15q
 - Prader-Willi syndrome: 15q
 - Cri-du chat syndrome: 5p
 - Wolf-Hirschhorn syndrome: 4p
 - Jacobsen syndrome: 11q
 - Langermeier syndrome: 8q
 - Hyperphagia
 - Obesity

National Institutes of Health. February 2014. <http://ghr.nlm.nih.gov>. Accessed February 11, 2015.




CONFIDENTIAL: Not for distribution, copyright reserved Quest Diagnostics 2015

Evolution of Screening:
MSS to NIPS

CONFIDENTIAL: Not for distribution, copyright reserved Quest Diagnostics 2015

Evolution of Fetal Aneuploidy Screening




Devers PL. J Genet Coun. 2013;22(3):291-295

* Nuchal translucency
** Applies to T21, T18, and T13
*** Some cases



CONFIDENTIAL: Not for distribution, copyright reserved Quest Diagnostics 2015

Current Options for Fetal Aneuploidy Screening



Discoff DA. Genet Med. 2008;11(11):818-821.



CONFIDENTIAL: Not for distribution, copyright reserved Quest Diagnostics 2015

Clinical Utility of NIPS: Action from Insight.

CONFIDENTIAL: Not for distribution, copyright reserved Quest Diagnostics 2015

Clinical Impact of NIPS

Total abnormal amniocentesis results = 2.3%

Aneuploidy; Robertsonian translocation 1.86%	
▪ T21	1.16%
▪ Sex chromosome aneuploidy	0.35%
▪ T18	0.23%
▪ T13	0.07%
▪ t(13:14)	0.05%

▪ Marker; unbalanced translocation	0.12%
▪ Other	0.09%

Balanced arrangements	0.20%
-----------------------	-------

97.7%
of amniocenteses have normal results

82% of total chromosomal abnormalities
Detectable now by NIPS*


9% of total chromosomal abnormalities
Reported detectable in literature by MPSS**

9% of total chromosomal abnormalities
Unlikely detectable by MPSS**

Ferguson-Smith MA. Prenat Diagn. 1984 Spring;4 Spec No:5-44.
CONFIDENTIAL: Not for distribution, copyright reserved Quest Diagnostics 2015



Clinical Impact of NIPS



McCullough et al. PLoS One, 2014
CONFIDENTIAL: Not for distribution, copyright reserved Quest Diagnostics 2015



Clinical Impact of NIPS on Invasive Prenatal Diagnosis

- Action from insight
 - Provides invaluable insight for informed decisions with regard to electing CVS or amniocentesis.
- Can be used as 1st tier screening (high-risk pregnancy) or 2nd tier screening (general and high-risk pregnancy)
- QNatal Advanced detection (sensitivity) rate = 100%
 - Validation series, confidence interval = 89.95% - 100%
 - General NIPS sensitivity rates range from 70% to 94%
 - Includes sex chromosome aneuploidy and microdeletions
- Higher positive predictive value (PPV) than MSS, NT, or a combination of both
- NIPS is screening; it is not diagnostic
 - False positives and negatives do occur; therefore, an invasive test is required to definitively confirm any positive result.

Benn PA. Obstet Gynecol. 2002;100(6):1168-1176.; Malone FD. N Eng J Med. 2005;353:2001-2011.; Swanson A. Curr Genet Med Rep. 2013;1(2):119-121.; Wang JC. Genet Med. 2014. doi:10.1038/gim.2014.92.
CONFIDENTIAL: Not for distribution, copyright reserved Quest Diagnostics 2015



Clinical Complications of NIPS: False Positives and Negatives


- Confined placental mosaicism (CPM)
- Co-twin demise/vanishing twin
- Fetal mosaicism
- Maternal mosaicism
- Maternal malignancy
- Collection/laboratory error
- Directly related to prevalence (positive and negative predictive value)

CONFIDENTIAL: Not for distribution, copyright reserved Quest Diagnostics 2015



Positive and Negative Predictive Value (PPV and NPV)

- Unlike sensitivity and specificity, PPV and NPV are dependent upon condition prevalence
 - PPV: percentage of patients with a positive test who truly have the condition
 - $PPV (\%) = \frac{\text{true positive}}{\text{true positive} + \text{false positive}} \times 100$
- PPV for NIPS is lower in the general pregnancy population, as the incidence of chromosomal aneuploidy is lower in these patients:



Wang JC. Genet Med. 2014. doi: 10.1038/gim.2014.92
CONFIDENTIAL: Not for distribution, copyright reserved Quest Diagnostics 2015



Summary of ACOG, ISPD, and NSGC Recommendations

- NIPS should be included in prenatal testing options offered to high-risk patient population
- NIPS should be offered within the context of informed consent, education, and counseling by a qualified provider
- NIPS positive results should be followed by recommending an option for confirmatory diagnostic testing
 - Negative results do not ensure an unaffected fetus

American College of Obstetricians and Gynecologists Committee on Genetics. Obstet Gynecol. 2012;120:1532-1534.
Devers PL. J Genet Coun. 2013;22(3):291-295.
Benn P. Prenat Diagn. 2012;32(1):1-2.
CONFIDENTIAL: Not for distribution, copyright reserved Quest Diagnostics 2015



QNatal Advanced
Performance Snapshot

CONFIDENTIAL: Not for distribution, copyright reserved Quest Diagnostics 2015



QNatal Advanced Test Content

Location	Syndrome Name	Estimated Frequency of Live Birth
Chromosome 21	Down	1/800
Chromosome 18	Edwards	1/6,000
Chromosome 13	Patau	1/7,000-33,000
Chromosome X	X, XXX, XXY, XYY	1/500
Sex Chromosomes		
22q*	DiGeorge	1/4,000
1p*	1p36	1/10,000
15q*	Angelman	1/20,000
15q*	Prader-Willi	1/20,000
5p*	Cri-du-chat	1/50,000
4p*	Wolf-Hirschhorn	1/50,000
11q*	Jacobsen	1/100,000
8q*	Langer-Giedion	Rare

*Microdeletion

National Institutes of Health. February 2014. <http://ghr.nlm.nih.gov>. Accessed February 11, 2015.

CONFIDENTIAL. Not for distribution. copyright reserved Quest Diagnostics 2015



Performance of QNatal Advanced

NIPS detects cell-free fetal DNA (cfDNA) of trophoblastic origin from the plasma of pregnant women


- 10% of DNA fragments in a pregnant woman's plasma are fetal

MPSS platform is used to detect maternal and fetal cfDNA

- MPSS can elucidate the chromosome of origin for a particular DNA fragment (see next slide)
- MPSS cannot differentiate between maternal and fetal DNA fragments

Bioinformatic analysis is used to determine the representation of fetal cfDNA

- GC normalization: routine correction to account for suboptimal amplification of GC-rich regions




Bianchi DW. *Nat Med*. 2012;18(7):1041-1051.; Bianchi DW et al. *Obstet Gynecol*. 2012;119(5):890-901.; Benn P. *Prenat Diagn*. 2012;32(1):1-2.; Zhao C et al. *PLoS One*. In press

CONFIDENTIAL. Not for distribution. copyright reserved Quest Diagnostics 2015




Performance of QNatal Advanced



CONFIDENTIAL. Not for distribution. copyright reserved Quest Diagnostics 2015



Performance of QNatal Advanced




CONFIDENTIAL. Not for distribution. copyright reserved Quest Diagnostics 2015



Validation of QNatal Advanced: Z-score Calculation

- Z-Score is used to determine positive/negative result for core trisomies
 - i.e., over-representation of chromosome material vs. appropriate representation of chromosome material
 - Z-Score = number of standard deviations from the mean (average)
 - Z-Score (chromosome 21) = $\frac{(\text{chromosome 21\%} - \text{mean chromosome 21\%})}{(\text{standard deviation chromosome 21\%})}$
- In assay development, all trisomies had Z-scores greater than 11 when fetal fraction was greater than 5%
- In validation, all unaffected pregnancies had Z-scores less than 4



Zhao C et al. PLoS One. In press
CONFIDENTIAL: Not for distribution, copyright reserved Quest Di.



Summary of QNatal Advanced

- Our laboratory-developed NIPS
- Can be performed as early as 10 weeks' gestation
 - At any gestational age thereafter without restriction
- Whole genome massively parallel shotgun sequencing (MPSS)
- Measures and reports fetal fraction as either sufficient or insufficient for analysis using a 5% cut-off
- Binary positive/negative result for core trisomies (T21, T18, T13)
- Clear reporting of sex chromosome aneuploidy and microdeletions as additional, incidental findings
- Opt out selections for microdeletions and fetal sex
- Validated for use in twin gestation
- 1-2% "no call" rate based on internal validation data

CONFIDENTIAL: Not for distribution, copyright reserved Quest Diagnostics 2015



QNatal Advanced

Sample Reports

CONFIDENTIAL: Not for distribution, copyright reserved Quest Diagnostics 2015

QNatal Advanced Sample Report: Clear and Concise Negative, Fetal Sex Opted Out Negative, Female

Report Header: Final EHSX, EHSXENG				Report Header: Final EHSX, EHSXENG			
Public Information		System Information		Client Information		Public Information	
TEST: EHSXENG	System: 00000001	Client P-175766	ADDRESS:	TEST: EHSXENG	System: 00000001	Client P-175766	ADDRESS:
Job: 00000076	AGE: 44	Prevalence: 001000	TEST CLIENT ID#:	Job: 00000076	AGE: 44	Prevalence: 001000	TEST CLIENT ID#:
Order: 1	Fetus: 19	Lab Ref: 00100000	TEST COLLECTOR ID#:	Order: 1	Fetus: 19	Lab Ref: 00100000	TEST COLLECTOR ID#:
Created: 07/13/2014 10:40:00Z	Received: 07/13/2014 10:40:00Z	Collected: 07/13/2014 10:40:00Z	COLLECTORVILLE, PA 19020-296	Collected: 07/13/2014 10:40:00Z	Received: 07/13/2014 10:40:00Z	Collected: 07/13/2014 10:40:00Z	COLLECTORVILLE, PA 19020-296
Received: 07/13/2014 10:40:00Z	Reported: 08/03/2014 10:28:00Z	Reported: 08/03/2014 10:28:00Z		Reported: 08/03/2014 10:28:00Z		Reported: 08/03/2014 10:28:00Z	
QNatal Advanced™				QNatal Advanced™			
Interpretation Summary				Interpretation Summary			
This specimen showed expected representation of chromosomes 21, 18, and 13 material.				This specimen showed expected representation of chromosomes 21, 18, 13 material.			
Chromosome Results		Fetal Sex Result		Chromosome Results		Fetal Sex Result	
Chromosome Tested	Results	Interpretation	Reference Range	Chromosome Tested	Results	Interpretation	Reference Range
Trisomy 21 (T21)	Negative	No apparent abnormality was detected. See "Limitation" below.	0-1	Trisomy 21 (T21)	Negative	No apparent abnormality was detected. See "Limitation" below.	0-1
Trisomy 18 (T18)	Negative	No apparent abnormality was detected. See "Limitation" below.	0-1	Trisomy 18 (T18)	Negative	No apparent abnormality was detected. See "Limitation" below.	0-1
Trisomy 13 (T13)	Negative	No apparent abnormality was detected. See "Limitation" below.	0-1	Trisomy 13 (T13)	Negative	No apparent abnormality was detected. See "Limitation" below.	0-1
Additional Chromosome Results		Interpretation		Additional Chromosome Results		Interpretation	
Sex Chromosomes	XX (Female)	No apparent abnormality was detected. See "Limitation" below.		Sex Chromosomes	XX (Female)	No apparent abnormality was detected. See "Limitation" below.	
Microdeletions	Not Detected	No apparent abnormality was detected. See "Limitation" below.		Microdeletions	Not Detected	No apparent abnormality was detected. See "Limitation" below.	
Laboratory Comments				Laboratory Comments			
See				See			
Client Provided High-Risk Indicators				Client Provided High-Risk Indicators			
Advanced Maternal Age Screening	Requested	Completed	Requested	Advanced Maternal Age Screening	Requested	Completed	Requested
Advanced Maternal Serum Screening	No	Preconception/Early History	No	Advanced Maternal Serum Screening	No	Preconception/Early History	No
<small> Results are preliminary unless stated post final Quest Diagnostics genetic consultation or call 800-521-6691 (ext. 4000) for assistance with interpretation of test results. </small>				<small> Results are preliminary unless stated post final Quest Diagnostics genetic consultation or call 800-521-6691 (ext. 4000) for assistance with interpretation of test results. </small>			

CONFIDENTIAL: Not for distribution, copyright reserved Quest Diagnostics 2015



QNatal Advanced Sample Report: Positive T21, Microdeletions Opted Out

Report Title: Fetal TEST, TESTENG

Patient Information	Specimen Information	Client Information
TEST, TESTENG DOB: 01/01/1978 Gender: F Race: NS Patient ID: 792	Specimen: 01/000010 Requestor: 001208 Lab Ref#: 00120901 Collected: 07/13/2014 10:18 EDT Received: 08/06/2014 10:18 EDT	Class #: 9702300 TEST CLIENT (DQ) 1301 COLLEGETVILLE RD COLLEGETVILLE, PA 19033-2998 Collected: 07/13/2014 10:18 EDT Received: 08/06/2014 10:18 EDT

QNatal Advanced™ Lab: 541

Interpretation Summary Lab: 541
This specimen showed additional representation of chromosome 21 material consistent with Trisomy 21 (Down syndrome). Further on-gene copy number and diagnostic testing is suggested. An expected representation of chromosomes 18 and 15 material was seen. Microdeletion testing was not performed per clinician request.

Chromosome Tested	Results	Findings	Concordance with a male fetus	Prevalence	Subtotal
Trisomy 21 (T21)	Positive	1 Chromosomal material detected	♂	1	1
Trisomy 18 (T18)	Negative	0 Chromosomal material detected		0	0
Trisomy 13 (T13)	Negative	0 Chromosomal material detected		0	0

Additional Chromosome Results

Chromosome Tested	Results	Interpretation
Sex Chromosomes	No anomaly	No apparent abnormality was detected. See "Limitations" below.
Microdeletions	Not performed	Microdeletion testing was not performed per clinician request.

Laboratory Comments
N/A

Clinician Provided High-Risk Indicators	Response	Component	Response
Advanced Maternal Age	Yes	Advanced (offered)	No
Advanced Maternal Serum Screening	No	Phenotypic/early history	No

Indicates possible pregnancies over true base frequency genetic carrier or cell loss (MPP) (see 10-103) for association with interpretation of these results.

QNatal Advanced Sample Report: XXY Male Finding

Report Title: Fetal TEST, TESTENG

Patient Information	Specimen Information	Client Information
TEST, TESTENG DOB: 01/01/1978 Gender: F Race: NS Patient ID: 792	Specimen: 01/000010 Requestor: 001208 Lab Ref#: 00120901 Collected: 07/13/2014 10:18 EDT Received: 08/06/2014 10:18 EDT	Class #: 9702300 TEST CLIENT (DQ) 1301 COLLEGETVILLE RD COLLEGETVILLE, PA 19033-2998 Collected: 07/13/2014 10:18 EDT Received: 08/06/2014 10:18 EDT

QNatal Advanced™ Lab: 541

Interpretation Summary Lab: 541
This specimen showed an abnormal representation of sex chromosome material (see "Sex Chromosomes" section below). Further on-gene copy number and diagnostic testing is suggested. An expected representation of chromosomes 21, 18, and 15 material was seen.

Chromosome Tested	Results	Findings	Concordance with a male fetus	Prevalence	Subtotal
Trisomy 21 (T21)	Negative	0 Chromosomal material detected	♂	1	1
Trisomy 18 (T18)	Negative	0 Chromosomal material detected		0	0
Trisomy 13 (T13)	Negative	0 Chromosomal material detected		0	0

Additional Chromosome Results

Chromosome Tested	Results	Interpretation
Sex Chromosomes	Trisomy XXY	Findings are suggestive of a 47,XXY chromosome karyotype, such as may be found in Klinefelter syndrome with associated symptoms.
Microdeletions	Not performed	No apparent abnormality was detected. See "Limitations" below.

Laboratory Comments
N/A

Clinician Provided High-Risk Indicators	Response	Component	Response
Advanced Maternal Age	Yes	Advanced (offered)	No
Advanced Maternal Serum Screening	No	Phenotypic/early history	No

Indicates possible pregnancies over true base frequency genetic carrier or cell loss (MPP) (see 10-103) for association with interpretation of these results.

CONFIDENTIAL - Not for distribution - copyright reserved Quest Diagnostics 2015



CONFIDENTIAL - Not for distribution - copyright reserved Quest Diagnostics 2015



QNatal Advanced Sample Report: Clear and Concise Negative, Male Twin(s) Negative, Female Twins

Report Title: Fetal TEST, TESTENG

Patient Information	Specimen Information	Client Information
TEST, TESTENG DOB: 01/01/1978 Gender: F Race: NS Patient ID: 792	Specimen: 01/000010 Requestor: 001208 Lab Ref#: 00120901 Collected: 07/13/2014 10:18 EDT Received: 08/06/2014 10:18 EDT	Class #: 9702300 TEST CLIENT (DQ) 1301 COLLEGETVILLE RD COLLEGETVILLE, PA 19033-2998 Collected: 07/13/2014 10:18 EDT Received: 08/06/2014 10:18 EDT

QNatal Advanced™ Lab: 541

Interpretation Summary Lab: 541
This specimen showed expected representation of chromosomes 21, 18, and 15 material.

Chromosome Tested	Results	Findings	Concordance with a male fetus	Prevalence	Subtotal
Trisomy 21 (T21)	Negative	0 Chromosomal material detected	♂	2	2
Trisomy 18 (T18)	Negative	0 Chromosomal material detected		0	0
Trisomy 13 (T13)	Negative	0 Chromosomal material detected		0	0

Additional Chromosome Results

Chromosome Tested	Results	Interpretation
Sex Chromosomes	No anomaly	No apparent abnormality was detected. See "Limitations" below.
Microdeletions	Not performed	No apparent abnormality was detected. See "Limitations" below.

Laboratory Comments
N/A

Clinician Provided High-Risk Indicators	Response	Component	Response
Advanced Maternal Age	Yes	Advanced (offered)	No
Advanced Maternal Serum Screening	No	Phenotypic/early history	No

Indicates possible pregnancies over true base frequency genetic carrier or cell loss (MPP) (see 10-103) for association with interpretation of these results.

CONFIDENTIAL - Not for distribution - copyright reserved Quest Diagnostics 2015



CONFIDENTIAL - Not for distribution - copyright reserved Quest Diagnostics 2015



Thank You QUESTIONS?